

Anti-Bag1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99933

Item Number: CM0002548

Introduction

Recombinant rabbit monoclonal antibody targeting human BAG1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Used in apoptosis and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	Bag1 Rabbit mAb
Remarks/Alias	HAP; BAG-1; RAP46; Bag1
Species	Human
Gene ID (Human)	573
Gene ID	573
Immunogen	Synthetic peptide, A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Bag1 (Q99933).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99933
Protein Weight	39kDa
Shipping	Ice bag