


Anti-Bag1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99933

Item Number: CM0015530



Introduction

Rabbit polyclonal antibody targeting human BAG1 (BAG-1/HAP/RAP46). Validated for WB, IF/ICC, and ELISA. Detects human and mouse BAG1 (~39 kDa). Ideal for apoptosis, chaperone, and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	BAG1 Rabbit pAb
Remarks/Alias	HAP; BAG-1; RAP46; BAG1
Species	Human
GeneID (Human)	573
GeneID	573
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99933
Protein Weight	39kDa
Shipping	Ice bag