

Anti-Bag3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 095817

Item Number: CM0014093

Introduction

BAG3 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets BAG family molecular chaperone regulator 3, a co-chaperone implicated in protein quality control and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	BAG3 Rabbit pAb
Remarks/Alias	BIS; MFM6; BAG-3; CAIR-1; BAG3
Species	Human
GeneID	9531
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 380-575 of human BAG3 (NP_004272.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095817
Protein Weight	62kDa
Shipping	Ice bag