

Anti-Bag4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O95429

Item Number: CM0019277



Introduction

Anti-BAG4 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets BAG4 (SODD), a molecular chaperone regulator inhibiting HSP70 and preventing constitutive TNFRSF1A signaling.

[Learn More](#)

Parameter	Value
Product Name	BAG4 Rabbit pAb
Remarks/Alias	SODD; BAG-4; BAG4
Species	Human
Gene ID (Human)	9530
Gene ID	9530
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human BAG4 (NP_001191807.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95429
Protein Weight	50kDa
Shipping	Ice bag