

Anti-Bag5 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9UI15

Item Number: CM0028062



Introduction

BAG5 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets the BAG5 co-chaperone of HSP70, essential for proteostasis and linked to cardiomyopathy. Ideal for cardiac and chaperone research.

[Learn More](#)

Parameter	Description
Product Name	BAG5 Rabbit pAb
Remarks/Alias	BAG-5; CMD2F; BAG5
Species	Human
Gene ID	9529
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 198-447 of human BAG5 (NP_004864.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UL15
Protein Weight	51kDa
Shipping	Ice bag