


Anti-Baiap2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uqb8

Item Number: CM0026780



Introduction

Rabbit polyclonal antibody targeting BAIAP2 (IRSp53/BAI-associated protein 2). Suitable for Western blot, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for actin cytoskeleton and neuronal development studies.

[Learn More](#)

Parameter	Value
Product Name	BAIAP2 Rabbit pAb
Remarks/Alias	BAP2, WAML, FLAF3, IRSP53, BAIAP2
Species	Human
Gene ID	10458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 362-521 of human BAIAP2 (NP_059344.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQB8
Protein Weight	61kDa
Shipping	Ice bag