


Anti-Baiap2L1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uhr4

Item Number: CM0027450



Introduction

Rabbit polyclonal antibody against BAIAP2L1 (IRTKS) validated for WB and ELISA. Detects human and mouse BAIAP2L1 protein, a BAR domain adapter involved in actin cytoskeleton reorganization. Suitable for cytoskeleton and infection research.

[Learn More](#)

Parameter	Value
Product Name	BAIAP2L1 Rabbit pAb
Remarks/Alias	IRTKS; BAIAP2L1
Species	Human
Gene ID (Human)	55971
Gene ID	55971
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 198-388 of human BAIAP2L1 (NP_061330.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHR4
Protein Weight	57kDa
Shipping	Ice bag