

Anti-Bak Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16611

Item Number: CM0006172



Introduction

Bak Rabbit mAb specifically detects human and monkey Bak. Validated for WB, IHC-P, and ELISA. A reliable choice for apoptosis research and IVD development, supported by CamelBio's one-stop raw material supply.

[Learn More](#)

Parameter	Value
Product Name	Bak Rabbit mAb
Remarks/Alias	BAK; CDN1; BCL2L7; BAK-LIKE; Bak
Species	Human
GeneID (Human)	578
GeneID	578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-82 of human Bak (NP_001179.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Monkey
SWISS	Q16611
Protein Weight	23kDa
Shipping	Ice bag