

Anti-Bak Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q16611

Item Number: CM0016726



Introduction

Rabbit polyclonal antibody targeting human Bak (BAK, BCL2L7) apoptosis regulator. Validated for WB, IF/ICC, IP, ELISA; reacts with human, mouse, rat. Ideal for mitochondrial apoptosis and Bcl-2 family research.

[Learn More](#)

Parameter	Value
Product Name	Bak Rabbit pAb
Remarks/Alias	BAK; CDN1; BCL2L7; BAK-LIKE; Bak
Species	Human
GeneID (Human)	578
GeneID	578
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human Bak (NP_001179.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16611
Protein Weight	23kDa
Shipping	Ice bag