

Anti-Banf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75531

Item Number: CM0005777



Introduction

High-quality rabbit monoclonal antibody against BANF1 (Barrier-to-autointegration factor), validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat BANF1. Suitable for chromatin organization, DNA damage, and viral integration research.

[Learn More](#)

Parameter	Value
Product Name	BANF1 Rabbit mAb
Remarks/Alias	BAF; NGPS; BCRP1; D14S1460; BANF1
Species	Human
GeneID (human)	8815
GeneID	8815
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-89 of human BANF1 (O75531)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75531
Protein Weight	10kDa
Shipping	Ice bag