

Anti-Banp Polyclonal Antibody For Wb, Elisa - Q8N9N5

Item Number: CM0030203



Introduction

Rabbit polyclonal antibody against human BANP (BEND1/SMAR1), validated for WB and ELISA. Detects human and mouse BANP, a chromatin regulator and tumor suppressor. Ideal for T-cell development, cell cycle, and splicing research.

[Learn More](#)

Parameter	Value
Product Name	BANP Rabbit pAb
Remarks/Alias	BEND1; SMAR1; SMARBP1; BANP
Species	Human
Gene ID (Human)	54971
Gene ID	54971
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human BANP (NP_060339.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N9N5
Protein Weight	56kDa
Shipping	Ice bag