


Anti-Bap1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q92560

Item Number: CM0026333



Introduction

BAP1 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat BAP1 (~80 kDa). Ideal for studying histone H2A deubiquitination and chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	BAP1 Rabbit pAb
Remarks/Alias	UBM2; UVM2; KURIS; TPDS1; UCHL2; hucep-6; HUCEP-13; BAP1
Species	Human
GeneID (Human)	8314
GeneID	8314
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human BAP1 (NP_004647.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92560
Protein Weight	80kDa
Shipping	Ice bag