

Anti-Bap31 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51572

Item Number: CM0006196



Introduction

BAP31 rabbit monoclonal antibody designed for WB, IHC-P, IF/ICC, ELISA. Detects BAP31 protein (28 kDa) in human, mouse, rat samples. Essential for studying ER function, protein transport, and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	BAP31 Rabbit mAb
Remarks/Alias	CDM; DDCH; BAP31; 6C6-AG; DELXQ28; DXS1357E; MICRODELXq28
Species	Human
Gene ID	10134
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-200 of human BAP31 (P51572)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51572
Protein Weight	28kDa
Shipping	Ice bag