

Anti-Bap31 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P51572

Item Number: CM0025517



Introduction

Rabbit polyclonal antibody against human BAP31 (BCAP31), an ER membrane protein involved in protein secretion and ERAD. Validated for WB, IF/ICC, IP, ELISA and cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	BAP31 Rabbit pAb
Remarks/Alias	CDM; DDCH; BAP31; 6C6-AG; DELXQ28; DXS1357E; MICRODELXq28
Species	Human
GeneID (Human)	10134
GeneID	10134
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 127-246 of human BAP31 (NP_001243376.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51572
Protein Weight	28kDa
Shipping	Ice bag