


Anti-Bat3/Bag6 Polyclonal Antibody For Wb, Elisa - P46379

Item Number: CM0020734



Introduction

Rabbit polyclonal antibody against human BAT3/BAG6, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studying protein quality control, apoptosis, and immune responses. UniProt P46379.

[Learn More](#)

Parameter	Value
Product Name	BAT3/BAG6 Rabbit pAb
Remarks/Alias	G3; BAT3; BAG-6; D6S52E; BAT3/BAG6
Species	Human
GeneID (human)	7917
GeneID	7917
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46379
Protein Weight	119kDa
Shipping	Ice bag