

Anti-Bat3/Bag6 Rabbit Monoclonal Antibody For Wb And Elisa - P46379

Item Number: CM0002961



Introduction

High-quality anti-BAT3/BAG6 rabbit monoclonal antibody validated for WB and ELISA. Detects human, mouse, rat BAG6, a chaperone involved in protein quality control and apoptosis. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	BAT3/BAG6 Rabbit mAb
Remarks/Alias	G3; BAT3; BAG-6; D6S52E; BAT3/BAG6
Species	Human
Gene ID (Human)	7917
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1033-1132 of human BAT3/BAG6 (P46379).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46379
Protein Weight	119kDa
Shipping	Ice bag