


Anti-Bax Polyclonal Antibody For Wb, If/Icc, Elisa - Q07812

Item Number: CM0013150



Introduction

Rabbit polyclonal antibody targeting human Bax (BCL2 Associated X). Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for apoptosis research. Immunogen: synthetic peptide (aa 1-100).

[Learn More](#)

Parameter	Value
Product Name	Bax Rabbit pAb
Remarks/Alias	BCL2 Associated X, Bcl-2-Like Protein 4, Bcl2-L-4, BCL2L4,BAX
Species	Human
Gene ID (Human)	581
Gene ID	581
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bax (NP_620116.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07812
Protein Weight	21kDa
Shipping	Ice bag