


Anti-Bax Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q07812

Item Number: CM0013232



Introduction

High-quality anti-Bax polyclonal antibody validated for Western blot, IHC-P, IF/ICC, and ELISA; reacts with human, mouse, and rat Bax. Ideal for apoptosis research.

[Learn More](#)

Parameter	Specification
Product Name	Bax Rabbit pAb
Remarks/Alias	BCL2 Associated X, Bcl-2-Like Protein 4, Bcl2-L-4, BCL2L4,BAX
Species	Human
GeneID (Human)	581
GeneID	581
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bax (NP_620116.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07812
Protein Weight	21kDa
Shipping	Ice bag