


Anti-Bax Rabbit Monoclonal Antibody (Ko Validated) For Wb, If/Icc, Ip, Elisa - Q07812

Item Number: CM0003437



Introduction

Validated Bax rabbit monoclonal antibody for apoptosis research. Suitable for WB, IF/ICC, IP, and ELISA in human, mouse, and rat. Targets SWISS Q07812 with a 21kDa protein weight.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Bax Rabbit mAb
Remarks/Alias	BCL2L4; ax
Species	Human
Gene ID (Human)	581
Gene ID	581
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07812
Protein Weight	21kDa
Shipping	Ice bag