


Anti-Bax Rabbit Polyclonal Antibody (Ko Validated) For Apoptosis Research - Q07812

Item Number: CM0012926



Introduction

KO validated rabbit polyclonal antibody against human Bax (BCL2 Associated X, BCL2L4). Detects 21 kDa protein. Suitable for WB, IHC-P, IF/ICC, IP, ELISA; cross-reacts with mouse, rat. Ideal for apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Bax Rabbit pAb
Remarks/Alias	BCL2 Associated X, Bcl-2-Like Protein 4, Bcl2-L-4, BCL2L4, BAX
Species	Human
GeneID (human)	581
GeneID	581
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bax (NP_620116.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07812
Protein Weight	21 kDa
Shipping	Ice bag