


Anti-Baz2A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uif9

Item Number: CM0022084



Introduction

High-quality rabbit polyclonal antibody targeting human BAZ2A (TIP5), validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for chromatin remodeling and rDNA silencing studies.

[Learn More](#)

Parameter	Value
Product Name	BAZ2A Rabbit pAb
Remarks/Alias	TIP5; WALp3; BAZ2A
Species	Human
GeneID (Human)	11176
GeneID	11176
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 97-295 of human BAZ2A (NP_038477.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIF9
Protein Weight	211kDa
Shipping	Ice bag