


Anti-Baz2A / Tip5 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uif9

Item Number: CM0021347



Introduction

High-quality rabbit polyclonal antibody against human BAZ2A / TIP5, validated for Western blot and ELISA. Specific to human and mouse targets, ideal for chromatin remodeling and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	BAZ2A / TIP5 Rabbit pAb
Remarks/Alias	TIP5; WALp3; BAZ2A / TIP5
Species	Human
GeneID (human)	11176
GeneID	11176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 97-295 of human BAZ2A (TIP5) (NP_038477.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UIF9
Protein Weight	211kDa
Shipping	Ice bag