

Anti-Bbs10 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Tam1

Item Number: CM0027511

Introduction

Anti-BBS10 rabbit polyclonal antibody (pAb) validated for Western blot (WB) and ELISA. Detects human and mouse BBS10 (~81 kDa). Ideal for studying ciliopathies, BBSome assembly, and molecular chaperone functions.

[Learn More](#)

Parameter	Value
Product Name	BBS10 Rabbit pAb
Remarks/Alias	C12orf58; BBS10
Species	Human
GeneID (Human)	79738
GeneID	79738
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 474-723 of human BBS10 (NP_078961.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TAM1
Protein Weight	81kDa
Shipping	Ice bag