

Anti-Bbs2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bxc9

Item Number: CM0026688

Introduction

Anti-BBS2 rabbit polyclonal antibody validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Useful for BBSome complex studies, ciliogenesis research, and Bardet-Biedl syndrome investigations.

[Learn More](#)

Parameter	Specification
Product Name	BBS2 Rabbit pAb
Remarks/Alias	BBS; RP74; BBS2
Species	Human
GeneID (Human)	583
GeneID	583
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-96 of human BBS2 (NP_114091.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BXC9
Protein Weight	80kDa
Shipping	Ice bag