

Anti-Bbs7 Rabbit Polyclonal Antibody For Wb/If/Icc/Elisa - Q8Iwz6

Item Number: CM0024634

Introduction

Rabbit polyclonal antibody targeting human BBS7 (BBS2L1), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Detects an ~80 kDa protein. Suitable for ciliopathy, BBSome complex, and SHH pathway research.

[Learn More](#)

Parameter	Value
Product Name	BBS7 Rabbit pAb
Remarks/Alias	BBS2L1; BBS7
Species	Human
GenelD (Human)	55212
GenelD	55212
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human BBS7 (NP_060660.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8Iwz6
Protein Weight	80kDa
Shipping	Ice bag