

Anti-Bc12 Polyclonal Antibody For Wb And Elisa - Q6Yul8

Item Number: CM0021922

Introduction

Rabbit polyclonal antibody targeting *Oryza sativa* BC12 (KIN4A), suitable for Western blot and ELISA. Validated for rice species. Useful for research on cell wall biogenesis and gibberellin signaling.

[Learn More](#)

Parameter	Value
Product Name	BC12 Rabbit pAb
Remarks/Alias	BC2; GDD1; KIF4-like; OsJ_28385; OJ1134_E08.39-1; OJ1134_E08.39-2; OJ1134_E08.39-3; BC12
Species	<i>Oryza sativa</i>
Gene ID	4346402
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 840-970 of <i>Oryza sativa</i> BC12. (Q6YUL8)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q6YUL8
Protein Weight	116kDa
Shipping	Ice bag