

## Anti-Bcar1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P56945

Item Number: CM0022482



### Introduction

BCAR1 Rabbit Polyclonal Antibody (P56945) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Targets p130cas, a central docking protein in cell adhesion and migration signaling pathways.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | BCAR1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | CAS; CAS1; CASS1; CRKAS; P130Cas; BCAR1                                                                            |
| Species          | Human                                                                                                              |
| Gene ID          | 9564                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 495-612 of human BCAR1 (NP_055382.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P56945                                                                                                             |
| Protein Weight   | 93kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |