

Anti-Bcas3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q9H6U6

Item Number: CM0028587



Introduction

Rabbit polyclonal antibody targeting human BCAS3, validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Suitable for angiogenesis, cancer, and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	BCAS3 Rabbit pAb
Remarks/Alias	MAAB, GAOB1, PHAF2, HEMARS, BCAS3
Species	Human
Gene ID (Human)	54828
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6U6
Protein Weight	101kDa
Shipping	Ice bag