


Anti-Bcat1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P54687

Item Number: CM0024216



Introduction

BCAT1 Rabbit pAb (P54687) is a polyclonal antibody validated for WB, IF/ICC, and ELISA. It reacts with human, mouse, and rat BCAT1, enabling studies on branched-chain amino acid metabolism and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	BCAT1 Rabbit pAb
Remarks/Alias	BCT1; PP18; BCATC; ECA39; MECA39; PNAS121; BCAT1
Species	Human
GeneID (Human)	586
GeneID	586
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-190 of human BCAT1 (NP_005495.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54687
Protein Weight	43kDa
Shipping	Ice bag