

Anti-Bcat2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15382

Item Number: CM0026724



Introduction

High-quality rabbit polyclonal antibody targeting human BCAT2 (mitochondrial branched-chain-amino-acid aminotransferase). Validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat. Ideal for metabolic research.

[Learn More](#)

Parameter	Value
Product Name	BCAT2 Rabbit pAb
Remarks/Alias	BCAM; BCT2; HVLI; PP18; BCATM; BCAT2
Species	Human
GeneID (human)	587
GeneID	587
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-392 of human BCAT2 (NP_001181.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15382
Protein Weight	44kDa
Shipping	Ice bag