

## Anti-Bccip Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9P287

Item Number: CM0028575



### Introduction

BCCIP rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat BCCIP. Ideal for DNA repair, cell cycle, and centrosome studies.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | BCCIP Rabbit pAb                                                                                                 |
| Remarks/Alias    | TOK1; TOK-1; BCCIP                                                                                               |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 56647                                                                                                            |
| GeneID           | 56647                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-314 of human BCCIP (NP_510868.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q9P287                                                                                                           |
| Protein Weight   | 36kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |