


Anti-Bcdin3D Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z5W3

Item Number: CM0021759



Introduction

High-quality rabbit polyclonal antibody targeting BCDIN3D (RNA 5'-monophosphate methyltransferase), validated for WB and ELISA. Cross-reacts with human and mouse. Suitable for studies on tRNA methylation and miRNA processing.

[Learn More](#)

Parameter	Value
Product Name	BCDIN3D Rabbit pAb
Remarks/Alias	BCDIN3D
Species	Human
GeneID (Human)	144233
GeneID	144233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human BCDIN3D (NP_859059.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z5W3
Protein Weight	33kDa
Shipping	Ice bag