


Anti-Bckdha Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P12694

Item Number: CM0030627



Introduction

Rabbit polyclonal antibody against human BCKDHA (MSU1, BCKDE1A), validated for WB, IF/ICC, ELISA. Detects human, mouse, rat. Targets mitochondrial branched-chain ketoacid dehydrogenase E1 alpha subunit, ~50kDa. Ideal for metabolism and maple syrup urine disease research.

[Learn More](#)

Parameter	Value
Product Name	BCKDHA Rabbit pAb
Remarks/Alias	MSU; MSUD1; OVD1A; BCKDE1A; BCKDHA
Species	Human
GeneID (Human)	593
GeneID	593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 276-445 of human BCKDHA (NP_000700.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12694
Protein Weight	50kDa
Shipping	Ice bag