


Anti-Bckdha Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P12694

Item Number: CM0018885



Introduction

Rabbit polyclonal antibody against human BCKDHA, validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying branched-chain amino acid metabolism and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	BCKDHA Rabbit pAb
Remarks/Alias	MSU; MSUD1; OVD1A; BCKDE1A; BCKDHA
Species	Human
Gene ID (Human)	593
Gene ID	593
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 276-445 of human BCKDHA (NP_000700.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12694
Protein Weight	50kDa
Shipping	Ice bag