


Anti-Bckdhb Monoclonal Antibody For Wb - P21953

Item Number: CM0006507



Introduction

Rabbit monoclonal antibody against BCKDHB (BCKDE1B), validated for WB and ELISA, reacts with human, mouse, rat. Used in research on branched-chain amino acid metabolism and maple syrup urine disease.

[Learn More](#)

Parameter	Value
Product Name	BCKDHB Rabbit mAb
Remarks/Alias	E1B; BCKDE1B; BCKDH E1-beta; BCKDHB
Species	Human
GeneID (Human)	594
GeneID	594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-392 of human BCKDHB (NP_898871.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21953
Protein Weight	43kDa
Shipping	Ice bag