


Anti-Bckdhb Rabbit Polyclonal Antibody For Wb, Elisa - P21953

Item Number: CM0012878



Introduction

Rabbit polyclonal antibody against human BCKDHB, validated for WB and ELISA, cross-reactive with mouse and rat. Targets mitochondrial 2-oxoisovalerate dehydrogenase subunit beta, involved in branched-chain amino acid catabolism.

[Learn More](#)

Parameter	Value
Product Name	BCKDHB Rabbit pAb
Remarks/Alias	E1B; BCKDE1B; BCKDH E1-beta; BCKDHB
Species	Human
GeneID (human)	594
GeneID	594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 183-392 of human BCKDHB (NP_898871.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P21953
Protein Weight	43kDa
Shipping	Ice bag