

Anti-Bckdk Monoclonal Antibody For Wb, Elisa - O14874

Item Number: CM0003687



Introduction

Rabbit monoclonal antibody targeting human BCKDK, validated for WB and ELISA. Cross-reacts with mouse and rat. Recognizes branched-chain alpha-ketoacid dehydrogenase kinase, a key regulator of BCAA and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	BCKDK Rabbit mAb
Remarks/Alias	BCK; BCKDKD; BCKDK
Species	Human
GeneID (human)	10295
GeneID	10295
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14874
Protein Weight	46kDa
Shipping	Ice bag