

Anti-Bckdk Polyclonal Antibody For Wb, If/Icc, Elisa - O14874

Item Number: CM0025693



Introduction

Rabbit polyclonal antibody targeting BCKDK (branched-chain alpha-ketoacid dehydrogenase kinase). Suitable for WB, IF/ICC, ELISA. Cross-reacts with human and rat. Ideal for research on BCAA metabolism, glucose and lipid regulation.

[Learn More](#)

Parameter	Value
Product Name	BCKDK Rabbit pAb
Remarks/Alias	BCK; BCKDKD; BCKDK
Species	Human
Gene ID (Human)	10295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 153-412 of human BCKDK (NP_005872.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O14874
Protein Weight	46kDa
Shipping	Ice bag