

Anti-Bcl-2 Monoclonal Antibody For Wb - P10415 / P49950

Item Number: CM0011669



Introduction

Anti-Bcl-2 rabbit monoclonal antibody validated for Western blot (WB), IHC-P, and ELISA. Recognizes human, mouse, and rat Bcl-2. Suitable for apoptosis, autophagy, and cancer research. Recombinant antibody for consistent performance.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Bcl-2 Rabbit PolymAb®
Remarks/Alias	Bcl-2; PPP1R50
Species	Human/Rat
GeneID (Human)	596
GeneID	596/24224
Immunogen	A synthetic peptide of human/rat Bcl2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10415/P49950
Protein Weight	26kDa
Shipping	Ice bag