

Anti-Bcl-2 Polyclonal Antibody For Wb, Ip, Elisa - P49950

Item Number: CM0022318

Introduction

Bcl-2 Rabbit pAb validated for WB, IP, ELISA with mouse cross-reactivity. Targets the N-terminal region of rat Bcl-2 (26 kDa), a key apoptosis regulator. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Bcl-2 Rabbit pAb
Remarks/Alias	Bcl-2,
Species	Rat
GeneID (Human)	596
GeneID	24224
Immunogen	Recombinant protein A synthetic peptide corresponding to a sequence within amino acids 1-100 of rat Bcl2(NP_058689.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse
SWISS	P49950
Protein Weight	26kDa
Shipping	Ice bag