

Anti-Bcl-2 Polyclonal Antibody For Wb - P10415

Item Number: CM0017298



Introduction

Rabbit polyclonal antibody against human Bcl-2, suitable for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Designed for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Bcl-2 Rabbit pAb
Remarks/Alias	Bcl-2; PPP1R50
Species	Human
GeneID (Human)	596
GeneID	596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of human Bcl-2 (NP_000624.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10415
Protein Weight	26kDa
Shipping	Ice bag