

Anti-Bcl-2 (Ppp1R50) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P10415

Item Number: CM0015400



Introduction

Rabbit polyclonal antibody against human Bcl-2 (PPP1R50), validated for WB, IHC-P, IF/ICC, IP, ELISA. Suppresses apoptosis; suitable for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Bcl-2 Rabbit pAb
Remarks/Alias	Bcl-2; PPP1R50
Species	Human
Gene ID (Human)	596
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bcl-2 (NP_000624.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10415
Protein Weight	26kDa
Shipping	Ice bag