

Anti-Bcl-2 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ihc-P, Ip, Elisa - P10415

Item Number: CM0003806



Introduction

High-quality Anti-Bcl-2 rabbit monoclonal antibody [KD Validated] validated for WB, IHC-P, IP, ELISA. Targets human Bcl-2 apoptosis regulator (P10415). Ideal for cancer research, apoptosis, and autophagy studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Bcl-2 Rabbit mAb
Remarks/Alias	Bcl-2; PPP1R50
Species	Human
GeneID (Human)	596
GeneID	596
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Bcl-2 (P10415)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	P10415
Protein Weight	26kDa
Shipping	Ice bag