


Anti-Bcl-W Polyclonal Antibody For Wb, Elisa - Q92843

Item Number: CM0018341



Introduction

Rabbit polyclonal antibody targeting human Bcl-W (BCL2L2), validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for apoptosis and cell survival research.

[Learn More](#)

Parameter	Value
Product Name	Bcl-W Rabbit pAb
Remarks/Alias	BCLW; BCL-W; PPP1R51; BCL2-L-2; Bcl-W
Species	Human
Gene ID (Human)	599
Gene ID	599
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-193 of human Bcl-W (NP_004041.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92843
Protein Weight	21kDa
Shipping	Ice bag