

Anti-Bcl-XL Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q07817

Item Number: CM0006436



Introduction

Rabbit monoclonal antibody against Bcl-XL, validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Ideal for apoptosis and mitochondrial research. Full product details, biological background, and sourcing support.

[Learn More](#)

Parameter	Specification
Product Name	Bcl-XL Rabbit mAb
Remarks/Alias	BCLX; BCL2L; Bcl-X; PPP1R52; BCL-XL/S; Bcl-XL
Species	Human
GeneID (Human)	598
GeneID	598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Bcl-XL (NP_612815.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07817
Protein Weight	26kDa
Shipping	Ice bag