

## Anti-Bcl-XL Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q07817

Item Number: CM0013433



### Introduction

Rabbit polyclonal antibody against Bcl-XL (aa 1-80). Validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for apoptosis and mitochondrial research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | Bcl-XL Rabbit pAb                                                                                                    |
| Remarks/Alias    | BCLX; BCL2L; Bcl-X; PPP1R52; BCL-XL/S; Bcl-XL                                                                        |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 598                                                                                                                  |
| Gene ID          | 598                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Bcl-XL (NP_001304848.1). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | Q07817                                                                                                               |
| Protein Weight   | 26 kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                              |