

Anti-Bcl-XL Rabbit Polyclonal Antibody (Kd Validated) For Wb, Elisa - Q07817

Item Number: CM0022801



Introduction

Rabbit polyclonal antibody against Bcl-XL, KD validated for Western blot and ELISA. Reacts with human, mouse, and rat Bcl-XL; ideal for apoptosis, cancer, and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] Bcl-XL Rabbit pAb
Remarks/Alias	BCLX; BCL2L; Bcl-X; PPP1R52; BCL-XL/S; [KD Validated] Bcl-XL
Species	Human
Gene ID (Human)	598
Gene ID	598
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Bcl-XL (NP_001304848.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07817
Protein Weight	26kDa
Shipping	Ice bag