


Anti-Bcl10 Monoclonal Antibody For Wb, Elisa - O95999

Item Number: CM0005924



Introduction

High-quality anti-Bcl10 rabbit monoclonal antibody validated for Western blot and ELISA. Specific for human Bcl10 (O95999), a key immune signaling adaptor. Ideal for NF-kappa-B and CBM complex research.

[Learn More](#)

Parameter	Value
Product Name	Bcl10 Rabbit mAb
Remarks/Alias	CLAP; mE10; CIPER; IMD37; c-E10; CARMEN; Bcl10
Species	Human
Gene ID (Human)	8915
Gene ID	8915
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95999
Protein Weight	26kDa
Shipping	Ice bag