


Anti-Bcl2A1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16548

Item Number: CM0017410



Introduction

BCL2A1 Rabbit Polyclonal Antibody, human-specific, validated for WB, IF/ICC, ELISA. Targets the hematopoietic anti-apoptotic protein Bfl-1 (BCL2A1) involved in cell survival and infection response. Ideal for apoptosis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	BCL2A1 Rabbit pAb
Remarks/Alias	GRS, ACC1, ACC2, BFL1, ACC-1, ACC-2, HBPA1, BCL2L5, BCL2A1
Species	Human
GeneID (Human)	597
GeneID	597
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 70-150 of human BCL2A1 (NP_004040.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q16548
Protein Weight	20kDa
Shipping	Ice bag