

Anti-Bcl2L12 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Hb09

Item Number: CM0003244

Introduction

BCL2L12 Rabbit mAb (Human, 37kDa, Q9HB09) validated for WB, IHC-P, IP, ELISA. Suitable for apoptosis research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	BCL2L12 Rabbit mAb
Species	Human
GeneID (Human)	83596
GeneID	83596
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q9HB09
Protein Weight	37kDa
Shipping	Ice bag